

LiveDemo

January 11, 2019

```
In [ ]: import warnings
```

```
from IPython.display import HTML
import matplotlib
from matplotlib import pyplot as plt
import numpy as np
plt.rcParams['figure.figsize'] = (20.0, 10.0)
matplotlib.rc('image', cmap='RdYlBu')

warnings.filterwarnings("ignore")
%matplotlib inline
```

```
In [ ]: #(Open dataset)
```

```
import karabo_data as kd
exmpl_file = '/gpfs/xfel/exp/XMPL/201750/p700000/proc/r0273/CORR-R0273-AGIPD03-S00000.h
hdf5_file = kd.H5File(exmpl_file)
```

Back to talk

```
In [ ]: #Live-demo (sel from train-id)
```

```
train_id, train_data = hdf5_file.train_from_id(198425246)
print(train_id, train_data.keys())
```

```
In [ ]: #Live-demo (sel from index)
```

```
train_id, train_data = hdf5_file.train_from_index(5)
print(train_id, train_data.keys())
```

```
In [ ]: #Live-demo (iteration)
```

```
for train_id, train_data in hdf5_file.trains():
    print(train_id, train_data.keys())
    break
```

Train data is of type dictionary. Hence the data can be accessed by giving keys:

```
In [ ]: #Live-demo show-dict
```

```
train_data['SPB_DET_AGIPD1M-1/DET/3CHO:xtdf']
```

```
In [ ]: #Live-demo data array and plot
```

```
print(train_data['SPB_DET_AGIPD1M-1/DET/3CHO:xtdf']['image.data'].shape)
```

```

print(type(train_data['SPB_DET_AGIPD1M-1/DET/3CHO:xtdf']['image.data']))
from matplotlib import pyplot as plt
fig = plt.figure(figsize=(10,10))
_ = plt.imshow(train_data['SPB_DET_AGIPD1M-1/DET/3CHO:xtdf']['image.data'][0].T,
                vmin=-50, vmax=100)

```

Back to talk

```

In [ ]: #Live-demo run dir
run_folder = '/gpfs/xfel/exp/XMPL/201750/p700000/raw/r0273'
run_dir = kd.RunDirectory(run_folder)
run_dir.info()

```

Back to talk

```

In [ ]: #Live-demo select with glob pattern
sel = run_dir.select('*/*DET/*', 'image.*')
sel.all_sources

```

Data can be access just like *previously* mentioned:

```

In [ ]: #Live-demo (iterator example)
for train_id, train_data in sel.trains():
    #print(train_id, train_data.keys())
    break
train_id, train_data = sel.train_from_index(5)
train_id, train_data.keys()

```

Back to talk

```

In [ ]: #Live-demo get series and plot
ph_flux = run_dir.get_series('SA1_XTD2_XGM/XGM/D00CS', 'pulseEnergy.photonFlux.value')
ph_flux.head()
fig = plt.figure(figsize=(10,5))
_ = ph_flux.plot()

```

Back to talk

```

In [ ]: #Live-demo get dataframe
fluxes_pos = run_dir.get_dataframe(fields=[("_XGM/*", "*.photonFlux"),
                                           ("*_XGM/*", "*.i[xy]Pos")])
fluxes_pos.head(10)

```

Back to talk

```

In [ ]: #Live-demo get array
xgm_intensity = run_dir.get_array('SA1_XTD2_XGM/XGM/D00CS:output', 'data.intensityTD',
                                  extra_dims=['pulseId'])
xgm_intensity

```

```
In [ ]: #Live-demo plot xgm data
fig = plt.figure(figsize=(10,8))
_= plt.imshow(xgm_intensity[:, :120].T,
               extent=(xgm_intensity.trainId[0], xgm_intensity.trainId[-1], 1, 120),
               origin='lower', cmap='RdYlBu_r')
```