

# OfflineAnalysis

January 23, 2019

```
In [2]: import warnings

        from IPython.display import HTML
        import matplotlib
        from matplotlib import pyplot as plt
        import numpy as np
        from matplotlib.cm import viridis
        import matplotlib
        viridis.set_bad('0.25',1)
        matplotlib.rc('image', cmap=viridis, origin='lower')

        warnings.filterwarnings("ignore")
        %matplotlib inline
```

A Closer Look at Detector Data:

Data from Modular Detectors is stored across multiple files

Using standard hdf5 tools can be tricky and tedious

## 0.1 *karabo-data*

- Python Library that supports Analysis of EuXFEL Data.
- It is Open Source and available on the maxwell cluster @ desy.

## 0.2 Aim of this presentation:

Give you a glimpse of what is possible and provide a ground base for your own data analysis at EuXFEL

## 0.3 Scenario I: Plot some detector data

An image from a single module of the AGIP-D Detector at SPB

- How could this data be retrieved and visualized?

```
In [ ]:
```

### 0.3.1 Reading a run directory:

One of the big advantages of karabo-data is that whole runs can be read with only *one* command:

```
In [3]: import karabo_data as kd
        run_dir = '/gpfs/xfel/exp/XMPL/201750/p700000/proc/r0273'
        run_data = kd.RunDirectory(run_dir)
```

Data can be accessed by: \* selecting trains by id's `.train_from_id` \* selecting trains by indexes `.train_from_index` \* iteration (looping) over trains `.trains`

**Let's select data based on indexes:**

```
In [4]: train_id, train_data = run_data.train_from_index(10)
        train_id, type(train_data)
```

```
Out[4]: (198425251, dict)
```

The data is stored in a so called *dictionary*. Hence the data can be accessed by giving keys:

```
In [5]: sorted(train_data.keys())
```

```
Out[5]: ['SPB_DET_AGIPD1M-1/DET/0CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/10CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/11CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/12CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/13CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/14CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/15CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/1CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/2CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/3CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/4CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/5CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/6CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/7CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/8CH0:xtdf',
         'SPB_DET_AGIPD1M-1/DET/9CH0:xtdf']
```

The properties of the above shown devices are can be accessed by selecting one of the names (keys):

```
In [6]: module_data = train_data['SPB_DET_AGIPD1M-1/DET/0CH0:xtdf']
        sorted(module_data.keys())
```

```
Out[6]: ['detector.data',
         'detector.trainId',
         'header.dataId',
         'header.linkId',
         'header.magicNumberBegin',
         'header.majorTrainFormatVersion',
```

```

'header.minorTrainFormatVersion',
'header.pulseCount',
'header.reserved',
'header.trainId',
'image.cellId',
'image.data',
'image.gain',
'image.length',
'image.mask',
'image.pulseId',
'image.status',
'image.trainId',
'metadata',
'trailer.checksum',
'trailer.magicNumberEnd',
'trailer.status',
'trailer.trainId']

```

The actual detector data is stored under the property value *image.data*

```

In [7]: img_ary = module_data['image.data']
        img_ary

```

```

Out[7]: array([[ 19.061886 ,  27.452316 ,  37.597824 , ...,  5.4725885 ,
                -7.8149734 , -25.363352 ],
               [ 33.70245 ,  22.299713 ,  0.93073153, ..., -4.7377877 ,
                -27.944302 , -5.8846965 ],
               [  2.7929337 , 31.081982 , 23.773945 , ..., -3.6882186 ,
                 1.8594275 , -9.089156 ],
               ...,
               [ -4.8543234 , -8.909287 , -8.339569 , ..., 24.457874 ,
                 43.014267 , 37.71819 ],
               [  6.3957534 , -1.3840029 ,  7.470259 , ..., 24.279312 ,
                 39.219994 , 47.97113 ],
               [ -5.3256435 , -8.172988 , 18.213696 , ..., 37.48085 ,
                 45.80402 , 87.84654 ]],

               [[185.36618 , 284.957 , 235.30476 , ..., 183.23979 ,
                 160.8983 , 121.14058 ],
               [274.3449 , 209.62767 , 173.04436 , ..., 141.82707 ,
                 63.43843 , 169.06215 ],
               [188.98767 , 250.20027 , 212.46991 , ..., 52.332073 ,
                 122.65631 , 133.87366 ],
               ...,
               [228.12357 , 198.33406 , 149.30414 , ..., 161.04172 ,
                 184.60376 , 191.79622 ],
               [295.8046 , 251.95976 , 82.540054 , ..., 358.52158 ,
                 212.06604 , 428.72366 ]],

```

[327.44974 , 183.78809 , 348.145 , ... , 239.44443 ,  
 223.869 , -13.2062025 ]],  
  
 [[350.2579 , 363.42407 , 374.40512 , ... , 254.54288 ,  
 383.2801 , 405.62738 ]],  
 [190.35628 , 185.2005 , 260.0179 , ... , 225.84837 ,  
 248.45126 , 335.21387 ]],  
 [181.68295 , 142.093 , 167.58374 , ... , 325.80933 ,  
 369.17328 , 389.54196 ]],  
 ... ,  
 [369.41043 , 337.76416 , 527.0745 , ... , 369.6595 ,  
 283.05963 , 485.79013 ]],  
 [555.3881 , 355.8697 , 286.35693 , ... , 543.9558 ,  
 389.25992 , 431.60147 ]],  
 [495.20718 , 418.86093 , 355.20557 , ... , 459.5037 ,  
 441.35126 , 443.27362 ]],  
  
 ... ,  
  
 [[ 12.515985 , 32.807766 , 17.720337 , ... , 116.88453 ,  
 120.97933 , 109.11631 ]],  
 [ 5.1309557 , 23.454872 , 10.613807 , ... , 87.08459 ,  
 87.28112 , 99.29501 ]],  
 [ 6.67852 , 8.346293 , 11.348902 , ... , 91.26027 ,  
 81.58455 , 105.83609 ]],  
 ... ,  
 [138.54237 , 125.0114 , 97.77815 , ... , 146.52628 ,  
 149.10728 , 167.16751 ]],  
 [129.00229 , 91.16923 , 92.66123 , ... , 161.34824 ,  
 153.23645 , 151.73201 ]],  
 [129.10573 , 111.126595 , 107.48337 , ... , 151.1517 ,  
 161.1087 , 317.03162 ]],  
  
 [[ 17.032906 , 23.841585 , 12.79358 , ... , 121.621025 ,  
 120.05122 , 112.59217 ]],  
 [ -1.0304956 , 9.197788 , 2.0586886 , ... , 82.65873 ,  
 101.2338 , 108.61823 ]],  
 [ -9.271917 , 1.0464283 , 13.519192 , ... , 92.814224 ,  
 112.34066 , 119.469475 ]],  
 ... ,  
 [131.70792 , 132.71964 , 119.65845 , ... , 153.29634 ,  
 166.48473 , 167.00223 ]],  
 [120.10571 , 113.27077 , 128.64536 , ... , 188.1266 ,  
 174.23666 , 175.4878 ]],  
 [123.17761 , 128.3278 , 136.94447 , ... , 167.48372 ,  
 163.49654 , 275.00546 ]],  
  
 [[159.35933 , 150.16685 , 156.23277 , ... , 264.5331 ,

```

276.90958 , 263.07755 ],
[167.7556 , 161.95264 , 161.65205 , ..., 249.0832 ,
254.59674 , 254.02673 ],
[163.67885 , 163.65152 , 164.29771 , ..., 246.56296 ,
251.98502 , 256.256 ],
...,
[309.0987 , 291.11078 , 274.75156 , ..., 318.2851 ,
309.00955 , 318.29727 ],
[283.04153 , 255.61102 , 297.74463 , ..., 325.1143 ,
314.2888 , 313.52383 ],
[308.6548 , 253.7317 , 272.84402 , ..., 318.43704 ,
311.13403 , 596.2051 ]]], dtype=float32)

```

The shape of the returned array is number of pulses x Y x X:

```
In [8]: img_ary.shape
```

```
Out[8]: (176, 512, 128)
```

To achieve the above shown plot three essential steps are necessary: \* stacking the detector data into one big Nd array \* loading the geometry description \* applying the loaded geometry to the data

Let's stack the data detector modules in the `module_data` dictionary first. The `stack_detector_data` method will create one array from all selected detector modules:

```
In [13]: train_img = kd.stack_detector_data(train_data, 'image.data', only='SPB_DET_AGIPD1M-1/DE
train_img.shape
```

```
Out[13]: (176, 16, 512, 128)
```

Now load the geometry information. This information is stored in a so called *geometry file* that describes the layout of the detector. Here we use the *.geom* format which is used by standard crystallography tools like CrystFEL. Karabo-data is able to handle this format:

```
In [14]: from karabo_data.geometry2 import AGIPD_1MGeometry
geom = AGIPD_1MGeometry.from_crystfel_geom('xfel.geom')
```


And let's apply the geometry and plot the data (11th pulse)

```
In [15]: img, center = geom.position_all_modules(train_img)
img.shape
```

```
Out[15]: (176, 1259, 1092)
```

```
In [16]: from matplotlib import pyplot as plt
plt.imshow(np.clip(img[10], -50, 1500))
```

```
Out[16]: <matplotlib.image.AxesImage at 0x2ad29c614a20>
```




Or even more easy you could *just* apply the `plot_data_fast` method:

```
In [17]: geom.plot_data_fast(np.clip(train_img[10], -50, 1500))
```

```
Out[17]:
```





### 0.3.2 Getting run information:

RunDirectory has a info method that displays a useful run experiment overview:

```
In [18]: run_folder = '/gpfs/exfel/exp/XMPL/201750/p700000/raw/r0273'
run_dir = kd.RunDirectory(run_folder)
run_dir.info()
```

```
# of trains:    156
Duration:      0:00:15.500000
```

First train ID: 198425241  
Last train ID: 198425396

16 detector modules (SPB\_DET\_AGIPD1M-1)  
e.g. module SPB\_DET\_AGIPD1M-1 0 : 512 x 128 pixels  
176 frames per train, 27456 total frames

2 instrument sources (excluding detectors):  
- SA1\_XTD2\_XGM/XGM/DOOCS:output  
- SPB\_XTD9\_XGM/XGM/DOOCS:output

13 control sources:  
- ACC\_SYS\_DOOCS/CTRL/BEAMCONDITIONS  
- SA1\_XTD2\_XGM/XGM/DOOCS  
- SPB\_IRU\_AGIPD1M/PSC/HV  
- SPB\_IRU\_AGIPD1M/TSENS/H1\_T\_EXTHOUS  
- SPB\_IRU\_AGIPD1M/TSENS/H2\_T\_EXTHOUS  
- SPB\_IRU\_AGIPD1M/TSENS/Q1\_T\_BLOCK  
- SPB\_IRU\_AGIPD1M/TSENS/Q2\_T\_BLOCK  
- SPB\_IRU\_AGIPD1M/TSENS/Q3\_T\_BLOCK  
- SPB\_IRU\_AGIPD1M/TSENS/Q4\_T\_BLOCK  
- SPB\_IRU\_AGIPD1M1/CTRL/MC1  
- SPB\_IRU\_AGIPD1M1/CTRL/MC2  
- SPB\_IRU\_VAC/GAUGE/GAUGE\_FR\_6  
- SPB\_XTD9\_XGM/XGM/DOOCS

Selecting data based by trains is simple with karabo-data but what if data should be selected across trains?

## 0.4 Scenario II: Extracting Data across trains with one value per train

<figcaption style="text-align: right">Photon flux <em>time-series</em> (by trainID)</figcaption>

- How can 1D data be extracted and plotted?

The `get_series` method can extract a series across trainID's for a given device and property:

The photon flux values are stored in the `pulseEnergy.photonFlux.value` property of device SA1\_XTD2\_XGM/XGM/DOOCS and has exactly one entry per train:

```
In [19]: ph_flux = run_dir.get_series('SA1_XTD2_XGM/XGM/DOOCS', 'pulseEnergy.photonFlux.value')
         type(ph_flux)
```

```
Out[19]: pandas.core.series.Series
```

```
In [20]: ph_flux.head(5)
```




```
Out[20]: trainId
198425241    500.519470
198425242    500.519470
198425243    502.727203
198425244    502.727203
198425245    504.070953
Name: SA1_XTD2_XGM/XGM/D00CS/pulseEnergy.photonFlux, dtype: float32
```

Pandas is a very useful data analysis library. More information is available under <https://pandas.pydata.org>.

```
In [21]: ph_flux.plot(figsize=(4,3))
```

```
Out[21]: <matplotlib.axes._subplots.AxesSubplot at 0x2ad29c6cf4e0>
```



#### 0.4.1 What if you wanted to get more than *one* device and/or property?

The `get_dataframe` method combines different sources into one single data object (also **pandas**):

```
In [22]: fluxes_pos = run_dir.get_dataframe(fields=[("*/XGM/D00CS", "*.i[xy]Pos")])
         type(fluxes_pos)
```

```
Out[22]: pandas.core.frame.DataFrame
```

```
In [23]: fluxes_pos.head(5)
```

```
Out[23]:          SPB_XTD9_XGM/XGM/D00CS/beamPosition.iyPos  \
trainId
198425241          -3.121433
```

198425242	-3.121433
198425243	-3.121433
198425244	-3.090523
198425245	-3.090523

	SPB_XTD9_XGM/XGM/D00CS/beamPosition.ixPos \
trainId	
198425241	5.512009
198425242	5.512009
198425243	5.512009
198425244	5.528512
198425245	5.528512

	SA1_XTD2_XGM/XGM/D00CS/beamPosition.iyPos \
trainId	
198425241	0.315761
198425242	0.315761
198425243	0.315761
198425244	0.341187
198425245	0.341187

	SA1_XTD2_XGM/XGM/D00CS/beamPosition.ixPos
trainId	
198425241	1.293711
198425242	1.293711
198425243	1.293711
198425244	1.336566
198425245	1.336566

## 0.5 Scenario III: Getting data with multiple values per train

X-ray gas intensity data is pulse resolved

- How can this 2D data be extracted and plotted?

the `get_array` method returns a data array that contains more than one value per train.

X-ray gas intensity data is pulse resolved and serves as an example:

the data can be accessed, similarly to `get_dataframe` by giving *device* name and *property*:

```
In [24]: xgm_intensity = run_dir.get_array('SA1_XTD2_XGM/XGM/D00CS:output', 'data.intensityTD')
xgm_intensity
```

```
Out[24]: <xarray.DataArray (trainId: 156, dim_0: 1000)>
array([[ 957.0532 , 1026.0005 ,  949.8755 , ...,  0.        ,  0.        ,
         0.        ],
       [ 763.8806 ,  794.2738 ,  868.2455 , ...,  0.        ,  0.        ,
         0.        ],
       [ 859.37   ,  995.1641 ,  838.5669 , ...,  0.        ,  0.        ,
         0.        ]],
```

```


...,
[ 945.2731 , 812.4336 , 839.45654, ..., 0. , 0. ,
 0. ],
[ 903.26855, 940.15125, 953.9436 , ..., 0. , 0. ,
 0. ],
[ 944.08386, 949.549 , 861.7509 , ..., 0. , 0. ,
 0. ]], dtype=float32)
Coordinates:
  * trainId (trainId) uint64 198425241 198425242 ... 198425395 198425396
Dimensions without coordinates: dim_0

In [25]: #Reset the default image size
matplotlib.rc('image', cmap='viridis', origin='lower')
matplotlib.rc('figure', figsize=(5,5))

In [26]: plt.imshow(xgm_intensity[:, :100].T)

Out[26]: <matplotlib.image.AxesImage at 0x2ad29c988eb8>

```



a labeled array (**xarray**) is returned. More information on labeled arrays can be found on <http://xarray.pydata.org>

## 0.6 Scenario IV: In the unlikely event if something goes wrong

A common source of errors is an invalid structure of the created data files. Hence a useful starting point to debug any errors is to check for valid file structures using the `karabo-data-validate` command

The command checks if:

- All .h5 files in a run can be opened, and the run contains at least one usable file.

- The list of train IDs in a file has no zeros except for padding at the end.
- Each train ID in a file is greater than the one before it.
- The indexes do not point to data beyond the end of a dataset.
- The indexes point to the start of the dataset, and then to successive chunks for successive trains, without gaps or overlaps between them.

## 1 How to install/get karabo-data?

Karabo-data is available on GitHub and there are multiple ways to install it: \* it is automatically available when you enter maxwell via jupyterhub <https://max-jhub.desy.de> \* it is installed in maxwell's anaconda environment module load anaconda/3 \* it can be install using *pip* `pip install (--user) karabo_data` \* the latest version could be downloaded from GitHub `git clone https://github.com/European-XFEL/karabo_data.git`

**A much more detailed documentation is available on *readthedocs* :**

<https://karabo-data.readthedocs.io/en/latest/>

In [ ]: