

```
In [2]: import warnings

from IPython.display import HTML
import matplotlib
from matplotlib import pyplot as plt
import numpy as np
from matplotlib.cm import viridis
import matplotlib
viridis.set_bad('0.25',1)
matplotlib.rc('image', cmap=viridis, origin='lower')

warnings.filterwarnings("ignore")
%matplotlib inline
```


OFFLINE DATA ANALYSIS

THE KARABO-DATA ECOSYSTEM

Martin Bergemann
Control and Analysis Software (CAS)




2019-01-24



```
ls /gpfs/exfel/exp/XMPL/201750/p700000/proc/r0273/
CORR-R0273-AGIPD00-S00000.h5  CORR-R0273-AGIPD05-S00002.h5  CORR-R0273-AGIPD10-S00004.h5
CORR-R0273-AGIPD00-S00001.h5  CORR-R0273-AGIPD05-S00003.h5  CORR-R0273-AGIPD10-S00005.h5
CORR-R0273-AGIPD00-S00002.h5  CORR-R0273-AGIPD05-S00004.h5  CORR-R0273-AGIPD11-S00000.h5
CORR-R0273-AGIPD00-S00003.h5  CORR-R0273-AGIPD05-S00005.h5  CORR-R0273-AGIPD11-S00001.h5
CORR-R0273-AGIPD00-S00004.h5  CORR-R0273-AGIPD06-S00000.h5  CORR-R0273-AGIPD11-S00002.h5
CORR-R0273-AGIPD00-S00005.h5  CORR-R0273-AGIPD06-S00001.h5  CORR-R0273-AGIPD11-S00003.h5
CORR-R0273-AGIPD01-S00000.h5  CORR-R0273-AGIPD06-S00002.h5  CORR-R0273-AGIPD11-S00004.h5
CORR-R0273-AGIPD01-S00001.h5  CORR-R0273-AGIPD06-S00003.h5  CORR-R0273-AGIPD11-S00005.h5
CORR-R0273-AGIPD01-S00002.h5  CORR-R0273-AGIPD06-S00004.h5  CORR-R0273-AGIPD12-S00000.h5
CORR-R0273-AGIPD01-S00003.h5  CORR-R0273-AGIPD06-S00005.h5  CORR-R0273-AGIPD12-S00001.h5
CORR-R0273-AGIPD01-S00004.h5  CORR-R0273-AGIPD07-S00000.h5  CORR-R0273-AGIPD12-S00002.h5
CORR-R0273-AGIPD01-S00005.h5  CORR-R0273-AGIPD07-S00001.h5  CORR-R0273-AGIPD12-S00003.h5
CORR-R0273-AGIPD02-S00000.h5  CORR-R0273-AGIPD07-S00002.h5  CORR-R0273-AGIPD12-S00004.h5
CORR-R0273-AGIPD02-S00001.h5  CORR-R0273-AGIPD07-S00003.h5  CORR-R0273-AGIPD12-S00005.h5
CORR-R0273-AGIPD02-S00002.h5  CORR-R0273-AGIPD07-S00004.h5  CORR-R0273-AGIPD13-S00000.h5
CORR-R0273-AGIPD02-S00003.h5  CORR-R0273-AGIPD07-S00005.h5  CORR-R0273-AGIPD13-S00001.h5
CORR-R0273-AGIPD02-S00004.h5  CORR-R0273-AGIPD08-S00000.h5  CORR-R0273-AGIPD13-S00002.h5
CORR-R0273-AGIPD02-S00005.h5  CORR-R0273-AGIPD08-S00001.h5  CORR-R0273-AGIPD13-S00003.h5
CORR-R0273-AGIPD03-S00000.h5  CORR-R0273-AGIPD08-S00002.h5  CORR-R0273-AGIPD13-S00004.h5
CORR-R0273-AGIPD03-S00001.h5  CORR-R0273-AGIPD08-S00003.h5  CORR-R0273-AGIPD13-S00005.h5
CORR-R0273-AGIPD03-S00002.h5  CORR-R0273-AGIPD08-S00004.h5  CORR-R0273-AGIPD14-S00000.h5
CORR-R0273-AGIPD03-S00003.h5  CORR-R0273-AGIPD08-S00005.h5  CORR-R0273-AGIPD14-S00001.h5
CORR-R0273-AGIPD03-S00004.h5  CORR-R0273-AGIPD09-S00000.h5  CORR-R0273-AGIPD14-S00002.h5
CORR-R0273-AGIPD03-S00005.h5  CORR-R0273-AGIPD09-S00001.h5  CORR-R0273-AGIPD14-S00003.h5
CORR-R0273-AGIPD04-S00000.h5  CORR-R0273-AGIPD09-S00002.h5  CORR-R0273-AGIPD14-S00004.h5
CORR-R0273-AGIPD04-S00001.h5  CORR-R0273-AGIPD09-S00003.h5  CORR-R0273-AGIPD14-S00005.h5
CORR-R0273-AGIPD04-S00002.h5  CORR-R0273-AGIPD09-S00004.h5  CORR-R0273-AGIPD15-S00000.h5
CORR-R0273-AGIPD04-S00003.h5  CORR-R0273-AGIPD09-S00005.h5  CORR-R0273-AGIPD15-S00001.h5
CORR-R0273-AGIPD04-S00004.h5  CORR-R0273-AGIPD10-S00000.h5  CORR-R0273-AGIPD15-S00002.h5
CORR-R0273-AGIPD04-S00005.h5  CORR-R0273-AGIPD10-S00001.h5  CORR-R0273-AGIPD15-S00003.h5
CORR-R0273-AGIPD05-S00000.h5  CORR-R0273-AGIPD10-S00002.h5  CORR-R0273-AGIPD15-S00004.h5
CORR-R0273-AGIPD05-S00001.h5  CORR-R0273-AGIPD10-S00003.h5  CORR-R0273-AGIPD15-S00005.h5
```

A Closer Look at Detector Data:



Data from Modular Detectors is stored across multiple files

```
In [4]: !lsxfel /gpfs/exfel/exp/XMPL/201750/p700000/raw/r0273
```

```
r0273 : Run directory
```

```
# of trains:      156
Duration:         0:00:15.500000
First train ID: 198425241
Last train ID: 198425396
```

```
16 detector modules (SPB_DET_AGIPD1M-1)
   e.g. module SPB_DET_AGIPD1M-1 0 : 512 x 128 pixels
   176 frames per train, 27456 total frames
```

```
2 instrument sources (excluding detectors):
- SA1_XTD2_XGM/XGM/D00CS:output
- SPB_XTD9_XGM/XGM/D00CS:output
```

```
13 control sources:
- ACC_SYS_D00CS/CTRL/BEAMCONDITIONS
- SA1_XTD2_XGM/XGM/D00CS
- SPB_IRU_AGIPD1M/PSC/HV
- SPB_IRU_AGIPD1M/TSSENS/H1_T_EXTHOUS
- SPB_IRU_AGIPD1M/TSSENS/H2_T_EXTHOUS
- SPB_IRU_AGIPD1M/TSSENS/Q1_T_BLOCK
- SPB_IRU_AGIPD1M/TSSENS/Q2_T_BLOCK
- SPB_IRU_AGIPD1M/TSSENS/Q3_T_BLOCK
- SPB_IRU_AGIPD1M/TSSENS/Q4_T_BLOCK
- SPB_IRU_AGIPD1M1/CTRL/MC1
- SPB_IRU_AGIPD1M1/CTRL/MC2
- SPB_IRU_VAC/GAUGE/GAUGE_FR_6
- SPB_XTD9_XGM/XGM/D00CS
```


Using standard hdf5 tools can be tricky and tedious

karabo-data

- Python Library that supports Analysis of EuXFEL Data.
- It is Open Source and available on the maxwell cluster @ desy.

Aim of this presentation: Give you a glimpse of what is possible and provide a ground base for your own data analysis at EuXFEL

Scenario I: Plot some detector data



An image from a single module of the AGIP-D Detector at SPB

- How could this data be retrieved and visualized?

Reading a run directory:

One of the big advantages of karabo-data is that whole runs can be read with only *one* command:

```
In [3]: import karabo_data as kd
run_dir = '/gpfs/exfel/exp/XMPL/201750/p700000/proc/r0273'
run_data = kd.RunDirectory(run_dir)
```

Data can be accessed by:

- selecting trains by id's → `.train_from_id`
- selecting trains by indexes → `.train_from_index`
- iteration (looping) over trains → `.trains`

Let's select data based on indexes:

```
In [4]: train_id, train_data = run_data.train_from_index(10)
train_id, type(train_data)
```

```
Out[4]: (198425251, dict)
```

The data is stored in a so called *dictionary*. Hence the data can be accessed by giving keys:

```
In [5]: sorted(train_data.keys())
```

```
Out[5]: ['SPB_DET_AGIPD1M-1/DET/0CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/10CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/11CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/12CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/13CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/14CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/15CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/1CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/2CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/3CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/4CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/5CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/6CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/7CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/8CH0:xtdf',  
         'SPB_DET_AGIPD1M-1/DET/9CH0:xtdf']
```

The properties of the above shown devices are can be accessed by selecting one of the names (keys):

```
In [6]: module_data = train_data['SPB_DET_AGIPD1M-1/DET/0CH0:xtdf']  
        sorted(module_data.keys())
```

```
Out[6]: ['detector.data',  
         'detector.trainId',  
         'header.dataId',  
         'header.linkId',  
         'header.magicNumberBegin',  
         'header.majorTrainFormatVersion',  
         'header.minorTrainFormatVersion',  
         'header.pulseCount',  
         'header.reserved',  
         'header.trainId',  
         'image.cellId',  
         'image.data',  
         'image.gain',  
         'image.length',  
         'image.mask',  
         'image.pulseId',  
         'image.status',  
         'image.trainId',  
         'metadata',  
         'trailer.checksum',  
         'trailer.magicNumberEnd',  
         'trailer.status',  
         'trailer.trainId']
```

The actual detector data is stored under the property value *image.data*

```
In [7]: img_ary = module_data['image.data']  
img_ary
```

```

Out[7]: array([[ 19.061886 ,  27.452316 ,  37.597824 , ...,  5.4725885 ,
                -7.8149734 , -25.363352 ],
               [ 33.70245 ,  22.299713 ,  0.93073153, ..., -4.7377877 ,
                -27.944302 , -5.8846965 ],
               [  2.7929337 ,  31.081982 ,  23.773945 , ..., -3.6882186 ,
                1.8594275 , -9.089156 ],
               ...,
               [ -4.8543234 , -8.909287 , -8.339569 , ...,  24.457874 ,
                43.014267 ,  37.71819 ],
               [  6.3957534 , -1.3840029 ,  7.470259 , ...,  24.279312 ,
                39.219994 ,  47.97113 ],
               [ -5.3256435 , -8.172988 ,  18.213696 , ...,  37.48085 ,
                45.80402 ,  87.84654 ]],

               [[185.36618 ,  284.957 ,  235.30476 , ...,  183.23979 ,
                160.8983 ,  121.14058 ],
               [274.3449 ,  209.62767 ,  173.04436 , ...,  141.82707 ,
                63.43843 ,  169.06215 ],
               [188.98767 ,  250.20027 ,  212.46991 , ...,  52.332073 ,
                122.65631 ,  133.87366 ],
               ...,
               [228.12357 ,  198.33406 ,  149.30414 , ...,  161.04172 ,
                184.60376 ,  191.79622 ],
               [295.8046 ,  251.95976 ,  82.540054 , ...,  358.52158 ,
                212.06604 ,  428.72366 ],
               [327.44974 ,  183.78809 ,  348.145 , ...,  239.44443 ,
                223.869 , -13.2062025 ]],

               [[350.2579 ,  363.42407 ,  374.40512 , ...,  254.54288 ,
                383.2801 ,  405.62738 ],
               [190.35628 ,  185.2005 ,  260.0179 , ...,  225.84837 ,
                248.45126 ,  335.21387 ],
               [181.68295 ,  142.093 ,  167.58374 , ...,  325.80933 ,
                369.17328 ,  389.54196 ],
               ...,
               [369.41043 ,  337.76416 ,  527.0745 , ...,  369.6595 ,
                283.05963 ,  485.79013 ],
               [555.3881 ,  355.8697 ,  286.35693 , ...,  543.9558 ,
                389.25992 ,  431.60147 ],
               [495.20718 ,  418.86093 ,  355.20557 , ...,  459.5037 ,
                441.35126 ,  443.27362 ]],

               ...,

               [[ 12.515985 ,  32.807766 ,  17.720337 , ...,  116.88453 ,
                120.97933 ,  109.11631 ],
               [  5.1309557 ,  23.454872 ,  10.613807 , ...,  87.08459 ,
                87.28112 ,  99.29501 ],
               [  6.67852 ,  8.346293 ,  11.348902 , ...,  91.26027 ,
                81.58455 ,  105.83609 ],
               ...,
               [138.54237 ,  125.0114 ,  97.77815 , ...,  146.52628 ,
                149.10728 ,  167.16751 ],
               [129.00229 ,  91.16923 ,  92.66123 , ...,  161.34824 ,
                153.23645 ,  151.73201 ],
               [129.10573 ,  111.126595 ,  107.48337 , ...,  151.1517 ,
                161.1087 ,  317.03162 ]],

               [[ 17.032906 ,  23.841585 ,  12.79358 , ...,  121.621025 ,
                120.05122 ,  112.59217 ],
               [ -1.0304956 ,  9.197788 ,  2.0586886 , ...,  82.65873 ,
                101.2338 ,  108.61823 ],
               [ -9.271917 ,  1.0464283 ,  13.519192 , ...,  92.814224 ,
                112.34066 ,  119.469475 ],
               ...,
               [131.70792 ,  132.71964 ,  119.65845 , ...,  153.29634 ,
                166.48473 ,  167.00223 ],
               [120.10571 ,  113.27077 ,  128.64536 , ...,  188.1266 ,
                174.23666 ,  175.4878 ],
               [123.17761 ,  128.3278 ,  136.94447 , ...,  167.48372 ,

```


The shape of the returned array is number of pulses x Y x X:

```
In [8]: img_ary.shape  
Out[8]: (176, 512, 128)
```

To achieve the above shown plot three essential steps are necessary:

- stacking the detector data into one big Nd array
- loading the geometry description
- applying the loaded geometry to the data

Let's stack the data detector modules in the `module_data` dictionary first. The `stack_detector_data` method will create one array from all selected detector modules:

```
In [13]: train_img = kd.stack_detector_data(train_data, 'image.data', only='SPB_DET_AGIPD1M-1/DET')  
train_img.shape  
Out[13]: (176, 16, 512, 128)
```

Now load the geometry information. This information is stored in a so called *geometry file* that describes the layout of the detector. Here we use the *.geom* format which is used by standard crystallography tools like CrystFEL. Karabo-data is able to handle this format:


```
In [14]: from karabo_data.geometry2 import AGIPD_1MGeometry  
geom = AGIPD_1MGeometry.from_crystfel_geom('xfel.geom')
```

And let's apply the geometry and plot the data (11th pulse)

```
In [15]: img, center = geom.position_all_modules(train_img)  
img.shape  
Out[15]: (176, 1259, 1092)
```

```
In [16]: from matplotlib import pyplot as plt  
plt.imshow(np.clip(img[10], -50, 1500))
```


```
Out[16]: <matplotlib.image.AxesImage at 0x2ad29c614a20>
```



Or even more easy you could *just* apply the `plot_data_fast` method:

```
In [17]: geom.plot_data_fast(np.clip(train_img[10], -50, 1500))
```

Out[17]:



Getting run information:

`RunDirectory` has a `info` method that displays a useful run experiment overview:

```
In [18]: run_folder = '/gpfs/exfel/exp/XMPL/201750/p700000/raw/r0273'
run_dir = kd.RunDirectory(run_folder)
run_dir.info()
```

```
# of trains:      156
Duration:         0:00:15.500000
First train ID: 198425241
Last train ID: 198425396


16 detector modules (SPB_DET_AGIPD1M-1)
  e.g. module SPB_DET_AGIPD1M-1 0 : 512 x 128 pixels
  176 frames per train, 27456 total frames

2 instrument sources (excluding detectors):
  - SA1_XTD2_XGM/XGM/D00CS:output
  - SPB_XTD9_XGM/XGM/D00CS:output

13 control sources:
  - ACC_SYS_D00CS/CTRL/BEAMCONDITIONS
  - SA1_XTD2_XGM/XGM/D00CS
  - SPB_IRU_AGIPD1M/PSC/HV
  - SPB_IRU_AGIPD1M/TSSENS/H1_T_EXTHOUS
  - SPB_IRU_AGIPD1M/TSSENS/H2_T_EXTHOUS
  - SPB_IRU_AGIPD1M/TSSENS/Q1_T_BLOCK
  - SPB_IRU_AGIPD1M/TSSENS/Q2_T_BLOCK
  - SPB_IRU_AGIPD1M/TSSENS/Q3_T_BLOCK
  - SPB_IRU_AGIPD1M/TSSENS/Q4_T_BLOCK
  - SPB_IRU_AGIPD1M1/CTRL/MC1
  - SPB_IRU_AGIPD1M1/CTRL/MC2
  - SPB_IRU_VAC/GAUGE/GAUGE_FR_6
  - SPB_XTD9_XGM/XGM/D00CS
```

Selecting data based by trains is simple with karabo-data but what if data should be selected across trains?

Scenario II: Extracting Data across trains with one value per train



Photon flux time-series (by trainID)

- How can 1D data be extracted and plotted?

The `get_series` method can extract a series across trainID's for a given device and property:

The photon flux values are stored in the `pulseEnergy.photonFlux.value` property of device `SA1_XTD2_XGM/XGM/D00CS` and has exactly one entry per train:

```
In [19]: ph_flux = run_dir.get_series('SA1_XTD2_XGM/XGM/D00CS', 'pulseEnergy.photonFlux.value')
         type(ph_flux)
```

```
Out[19]: pandas.core.series.Series
```


```
In [20]: ph_flux.head(5)
```

```
Out[20]: trainId
198425241    500.519470
198425242    500.519470
198425243    502.727203
198425244    502.727203
198425245    504.070953
Name: SA1_XTD2_XGM/XGM/D00CS/pulseEnergy.photonFlux, dtype: float32
```

Pandas is a very useful data analysis library. More information is available under <https://pandas.pydata.org> (<https://pandas.pydata.org>).

```
In [21]: ph_flux.plot(figsize=(4,3))
```

```
Out[21]: <matplotlib.axes._subplots.AxesSubplot at 0x2ad29c6cf4e0>
```



What if you wanted to get more than *one* device and/or property?

The `get_dataframe` method combines different sources into one single data object (also **pandas**):

```
In [22]: fluxes_pos = run_dir.get_dataframe(fields=[("*/XGM/D00CS", "*.i[xy]Pos")
         ])
         type(fluxes_pos)
```

```
Out[22]: pandas.core.frame.DataFrame
```



In [23]:

fluxes_pos.head(5)

Out[23]:

	SPB_XTD9_XGM/XGM /DOOCS /beamPosition.iyPos	SPB_XTD9_XGM/XGM /DOOCS /beamPosition.ixPos	SA1_XTD2_XGM/XGM /DOOCS /beamPosition.iyPos	SA1_XTD2_XGM/XGM /DOOCS /beamPosition.ixPos
trainId				
198425241	-3.121433	5.512009	0.315761	1.293711
198425242	-3.121433	5.512009	0.315761	1.293711
198425243	-3.121433	5.512009	0.315761	1.293711
198425244	-3.090523	5.528512	0.341187	1.336566
198425245	-3.090523	5.528512	0.341187	1.336566

Scenario III: Getting data with multiple values per train



- How can this 2D data be extracted and plotted?

the `get_array` method returns a data array that contains more than one value per train.

X-ray gas intensity data is pulse resolved and serves as an example:

the data can be accessed, similarly to `get_dataframe` by giving *device* name and *property*:


```
In [24]: xgm_intensity = run_dir.get_array('SA1_XTD2_XGM/XGM/D00CS:output', 'data
.intensityTD')
xgm_intensity
```

```
Out[24]: <xarray.DataArray (trainId: 156, dim_0: 1000)>
array([[ 957.0532 , 1026.0005 ,  949.8755 , ...,  0.      ,  0.      ,
         0.      ],
       [ 763.8806 ,  794.2738 ,  868.2455 , ...,  0.      ,  0.      ,
         0.      ],
       [ 859.37   ,  995.1641 ,  838.5669 , ...,  0.      ,  0.      ,
         0.      ],
       ...,
       [ 945.2731 ,  812.4336 ,  839.45654, ...,  0.      ,  0.      ,
         0.      ],
       [ 903.26855,  940.15125,  953.9436 , ...,  0.      ,  0.      ,
         0.      ],
       [ 944.08386,  949.549  ,  861.7509 , ...,  0.      ,  0.      ,
         0.      ]], dtype=float32)
Coordinates:
  * trainId  (trainId) uint64 198425241 198425242 ... 198425395 198425396
Dimensions without coordinates: dim_0
```

```
In [25]: #Reset the default image size
matplotlib.rc('image', cmap='viridis', origin='lower')
matplotlib.rc('figure', figsize=(5,5))
```

```
In [26]: plt.imshow(xgm_intensity[:, :100].T)
```

```
Out[26]: <matplotlib.image.AxesImage at 0x2ad29c988eb8>
```



a labeled array (**xarray**) is returned. More information on labeled arrays can be found on <http://xarray.pydata.org> (<http://xarray.pydata.org>)

Scenario IV: In the unlikely event if something goes wrong

A common source of errors is an invalid structure of the created data files. Hence a useful starting point to debug any errors is to check for valid file structures using the `karabo-data-validate` command

```
$: karabo-data-validate /gpfs/exfel/exp/XMPL/201750/p700000/raw/r0273
Checking run directory: /gpfs/exfel/exp/XMPL/201750/p700000/raw/r0273
No problems found
```

The command checks if:

- All .h5 files in a run can be opened, and the run contains at least one usable file.
- The list of train IDs in a file has no zeros except for padding at the end.
- Each train ID in a file is greater than the one before it.
- The indexes do not point to data beyond the end of a dataset.
- The indexes point to the start of the dataset, and then to successive chunks for successive trains, without gaps or overlaps between them.

How to install/get karabo-data?

Karabo-data is available on GitHub and there are multiple ways to install it:

- it is automatically available when you enter maxwell via jupyterhub → <https://max-jhub.desy.de> (<https://max-jhub.desy.de>)
- it is installed in maxwell's anaconda environment → `module load anaconda/3`
- it can be install using `pip` → `pip install (--user) karabo_data`
- the latest version could be downloaded from GitHub → `git clone https://github.com/European-XFEL/karabo_data.git`

A much more detailed documentation is available on *readthedocs* :

<https://karabo-data.readthedocs.io/en/latest/> (<https://karabo-data.readthedocs.io/en/latest/>)

In []: