

KaraboData_Talk

January 11, 2019

```
In [1]: import warnings

        from IPython.display import HTML
        import matplotlib
        from matplotlib import pyplot as plt
        import numpy as np
        plt.rcParams['figure.figsize'] = (20.0, 10.0)
        matplotlib.rc('image', cmap='RdYlBu')

        warnings.filterwarnings("ignore")
        %matplotlib inline
```

1 Offline data analysis

```
<div class="intro-body">
  <div class="intro_h1"><h3>The karabo-data ecosystem</h3></div>
  <p><strong><span class="a">Martin Bergemann</span></strong><span class="b"></span><span>CAS
```

A Closer Look at Detector Data:

Data from Modular Detectors is stored across multiple files

Using standard hdf5 tools can be tricky and tedious

1.1 Solution: *karabo-data*

- Python Library that supports Analysis of EuXFEL Data.
- It is Open Source and available on the max-well cluster.

1.2 Aim:

Give you a glimpse of what is possible and provide a ground base for your own data analysis at EuXFEL

1.3 Scenario I: Take a look at some detector data

An image from a single module of the AGIP-D Detector at SPB

- How could this data be retrieved and visualized?

1.3.1 Reading single files (demo)

```
In [2]: #Live-demo
import karabo_data as kd
exmpl_file = '/gpfs/xfel/exp/XMPL/201750/p700000/proc/r0273/CORR-R0273-AGIPD03-S00000.h
hdf5_file = kd.H5File(exmpl_file)
```

Data can be accessed by: * selecting trains id's * selecting train indexes * iteration over trains
Let's select data based on trains: (demo)

```
In [3]: #Live-demo (sel from train-id)
train_id, train_data = hdf5_file.train_from_id(198425246)
print(train_id, train_data.keys())

198425246 dict_keys(['SPB_DET_AGIPD1M-1/DET/3CH0:xtdf'])
```

```
In [4]: #Live-demo (sel from index)
train_id, train_data = hdf5_file.train_from_index(5)
print(train_id, train_data.keys())

198425246 dict_keys(['SPB_DET_AGIPD1M-1/DET/3CH0:xtdf'])
```

```
In [5]: #Live-demo (iteration)
for train_id, train_data in hdf5_file.trains():
    print(train_id, train_data.keys())
    break

198425241 dict_keys(['SPB_DET_AGIPD1M-1/DET/3CH0:xtdf'])
```

Train data is of type *dictionary*. Hence the data can be accessed by giving keys:

```
In [6]: train_data['SPB_DET_AGIPD1M-1/DET/3CH0:xtdf']

Out[6]: {'detector.data': array([ 0, 64, 114, ..., 0, 0, 0], dtype=uint8),
'detector.trainId': 198425241,
'header.dataId': 0,
'header.linkId': 18446744069414584335,
'header.magicNumberBegin': array([-50, -6, -17, -66, 70, 68, 84, 88], dtype=int8),
'header.majorTrainFormatVersion': 1,
'header.minorTrainFormatVersion': 0,
'header.pulseCount': 176,
'header.reserved': array([0, 0, 0, 0, 0, 0, 0, 0, 0, 0, 0, 0, 0, 0, 0, 0], dtype=uint8),
'header.trainId': 198425241,
'image.cellId': array([ 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12,
13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25,
26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38,
39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 49, 50, 51,
```

```

52, 53, 54, 55, 56, 57, 58, 59, 60, 61, 62, 63, 64,
65, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77,
78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90,
91, 92, 93, 94, 95, 96, 97, 98, 99, 100, 101, 102, 103,
104, 105, 106, 107, 108, 109, 110, 111, 112, 113, 114, 115, 116,
117, 118, 119, 120, 121, 122, 123, 124, 125, 126, 127, 128, 129,
130, 131, 132, 133, 134, 135, 136, 137, 138, 139, 140, 141, 142,
143, 144, 145, 146, 147, 148, 149, 150, 151, 152, 153, 154, 155,
156, 157, 158, 159, 160, 161, 162, 163, 164, 165, 166, 167, 168,
169, 170, 171, 172, 173, 174, 175], dtype=uint16),
'image.data': array([[ 2.86249847e+01,  2.84662476e+01,  2.94381580e+01, ...,
    3.97410750e+00,  1.53232050e+01,  2.65694761e+00],
 [ 3.76776047e+01,  1.66408386e+01,  1.25324402e+01, ...,
    3.53553891e+00, -7.76110840e+00,  5.06602287e+00],
 [ 2.33586578e+01,  1.33027906e+01,  1.29060669e+01, ...,
   -3.58545852e+00, -8.78809452e+00, -1.76954997e+00],
 ...,
 [ 2.67929325e+01,  3.82321243e+01,  1.69615307e+01, ...,
   -5.36543350e+01,  1.77964573e+01,  2.29760933e+01],
 [ 5.43048973e+01,  2.40606403e+01,  2.68685341e+01, ...,
    1.56021561e+02,  0.00000000e+00,  1.39021225e+02],
 [ 5.41300812e+01,  3.99643478e+01,  3.18235302e+01, ...,
    4.46845245e+00, -3.83902435e+02,  1.09092340e+01]],

[[ 2.82328766e+02,  2.13707367e+02,  3.85529510e+02, ...,
    2.87144348e+02,  2.49376450e+02,  3.46322205e+02],
 [ 1.72267563e+02,  2.91269714e+02,  2.23787582e+02, ...,
    1.87312103e+02,  1.91055298e+02,  1.92247986e+02],
 [ 1.74026535e+02,  2.09151901e+02,  2.23993912e+02, ...,
    1.28322266e+02,  6.25818977e+01,  1.35687103e+02],
 ...,
 [ 9.67379822e+02,  4.64296753e+02,  1.03795129e+03, ...,
    4.12897301e+01,  2.44968739e+01,  1.36369019e+02],
 [ 1.12769836e+03,  8.79650452e+02,  6.07432495e+02, ...,
    3.33366852e+02,  2.61182800e+02,  4.13849525e+01],
 [ 1.37194507e+03,  1.16530359e+03,  9.33655640e+02, ...,
   -1.40187941e+01, -1.80979218e+02, -1.05465247e+03]],

[[ 2.74763580e+02,  2.57296173e+02,  1.92696045e+02, ...,
    5.14890625e+02,  4.83632965e+02,  4.13015289e+02],
 [ 1.88023987e+02,  1.00196587e+02,  1.05712891e+02, ...,
    3.33557678e+02,  4.17891815e+02,  4.55165558e+02],
 [ 2.97341156e+02,  1.76079010e+02,  8.91075516e+01, ...,
    2.00919449e+02,  2.54712372e+02,  2.78320801e+02],
 ...,
 [ 1.41089551e+03,  1.34565967e+03,  1.08119263e+03, ...,
    7.16647387e+00,  2.45052700e+01,  5.89854355e+01],
 [ 1.80306104e+03,  1.42946960e+03,  9.75338379e+02, ...,

```

```

        1.60704651e+02, 2.65418762e+02, 1.36963837e+02],
    [ 2.03517444e+03, 1.21541614e+03, 9.28263855e+02, ...,
      9.67607212e+00, -2.76129272e+02, -5.88472351e+02]],

    ...,

    [[ 3.09872532e+00, -5.08774281e-01, 1.51738377e+01, ...,
        9.67595291e+01, 1.02846840e+02, 9.79162140e+01],
      [-2.02406406e+00, -8.11292076e+00, -1.02228737e+01, ...,
        1.00479095e+02, 1.04622551e+02, 1.17755341e+02],
      [ 1.53183994e+01, 9.43448830e+00, -9.73048973e+00, ...,
        9.80847092e+01, 9.90229797e+01, 9.15223618e+01],
      ...,
      [ 9.65496826e+01, 9.00805054e+01, 8.23390045e+01, ...,
        3.99469910e+01, -4.39730103e+02, 8.46993359e+03],
      [ 1.09837486e+02, 9.43753433e+01, 1.01021187e+02, ...,
        -7.42510234e+04, -1.19970758e+05, -7.43315469e+04],
      [ 1.11256569e+02, 9.82772827e+01, 1.06616699e+02, ...,
        7.79505127e+02, -1.85362646e+03, -1.54386367e+04]],

    ...,

    [[ 2.27756405e+01, 2.25000038e+01, 2.32766819e+01, ...,
        1.05306572e+02, 1.19960114e+02, 1.06827599e+02],
      [ 1.22687225e+01, 1.11412373e+01, -3.10237098e+00, ...,
        1.11624672e+02, 9.99300003e+01, 1.20338463e+02],
      [ 2.03352380e+00, 3.13974047e+00, 1.52829742e+01, ...,
        1.17073769e+02, 1.08744743e+02, 1.05334351e+02],
      ...,
      [ 1.06949814e+02, 9.66776505e+01, 9.08718491e+01, ...,
        -1.01162968e+01, -8.50073535e+03, -1.30835596e+04],
      [ 1.24415733e+02, 9.19042053e+01, 9.73192291e+01, ...,
        -2.65518613e+04, 1.87262070e+04, -1.64417363e+04],
      [ 1.32533249e+02, 1.10008453e+02, 1.15432304e+02, ...,
        3.16230347e+02, -8.02921484e+03, -1.48106921e+03]],

    ...,

    [[ 1.53365967e+02, 1.39493073e+02, 1.64242966e+02, ...,
        2.43825836e+02, 2.61614197e+02, 2.44643021e+02],
      [ 1.56373383e+02, 1.53281754e+02, 1.43108688e+02, ...,
        2.61418030e+02, 2.48436081e+02, 2.61712524e+02],
      [ 1.51726349e+02, 1.45641235e+02, 1.53762711e+02, ...,
        2.48970200e+02, 2.41733551e+02, 2.53827728e+02],
      ...,
      [ 2.72267700e+02, 2.46381104e+02, 2.55895477e+02, ...,
        1.01002051e+03, -2.68937812e+04, -6.32623047e+04],
      [ 2.82559753e+02, 2.68116241e+02, 2.57689056e+02, ...,
        -6.13691211e+03, -5.68883984e+03, -3.82133008e+03],
      [ 2.84443420e+02, 2.49587311e+02, 2.78960602e+02, ...,
        3.70983740e+03, 9.68601562e+03, 1.14233716e+03]]],
dtype=float32),

```

```

'image.gain': array([[0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    ...,
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 1, 0],
                    [0, 0, 0, ..., 0, 0, 0]],

                   [[0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    ...,
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 1]],

                   [[0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    ...,
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 1]],

                   ...,

                   [[0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    ...,
                    [0, 0, 0, ..., 1, 1, 1],
                    [0, 0, 0, ..., 2, 2, 2],
                    [0, 0, 0, ..., 1, 1, 2]],


                   [[0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    ...,
                    [0, 0, 0, ..., 1, 2, 2],
                    [0, 0, 0, ..., 2, 2, 2],
                    [0, 0, 0, ..., 1, 1, 2]],

                   [[0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    [0, 0, 0, ..., 0, 0, 0],
                    ...,
                    [0, 0, 0, ..., 1, 2, 2],
                    [0, 0, 0, ..., 2, 2, 2],

```



```
fig = plt.figure(figsize=(10,10))
_ = plt.imshow(train_data['SPB_DET_AGIPD1M-1/DET/3CH0:xtdf']['image.data'][0].T,
               vmin=-50, vmax=100)
```



1.4 Reading whole Runs:

The main advantage of karabo-data is that not only single files at a time but whole runs can be read:

```
In [13]: #Live-demo
run_folder = '/gpfs/xfel/exp/XMPL/201750/p700000/raw/r0273'
run_dir = kd.RunDirectory(run_folder)
run_dir.info()
```

```
# of trains:      156
Duration:         0:00:15.500000
First train ID: 198425241
Last train ID: 198425396
```

```
16 detector modules (SPB_DET_AGIPD1M-1)
  e.g. module SPB_DET_AGIPD1M-1 0 : 512 x 128 pixels
  176 frames per train, 27456 total frames
```

2 instrument sources (excluding detectors):

- SA1_XTD2_XGM/XGM/DOOCS:output
- SPB_XTD9_XGM/XGM/DOOCS:output

13 control sources:

- ACC_SYS_DOOCS/CTRL/BEAMCONDITIONS
- SA1_XTD2_XGM/XGM/DOOCS
- SPB_IRU_AGIPD1M/PSC/HV
- SPB_IRU_AGIPD1M/TSENS/H1_T_EXTHOUS
- SPB_IRU_AGIPD1M/TSENS/H2_T_EXTHOUS
- SPB_IRU_AGIPD1M/TSENS/Q1_T_BLOCK
- SPB_IRU_AGIPD1M/TSENS/Q2_T_BLOCK

- SPB_IRU_AGIPD1M/TSENS/Q3_T_BLOCK
- SPB_IRU_AGIPD1M/TSENS/Q4_T_BLOCK
- SPB_IRU_AGIPD1M1/CTRL/MC1
- SPB_IRU_AGIPD1M1/CTRL/MC2
- SPB_IRU_VAC/GAUGE/GAUGE_FR_6
- SPB_XTD9_XGM/XGM/DOOCS

1.5 Data subsets

Yet another powerful tool is the *select* method. It selects a subset of sources and keys from the run directory:

```
In [14]: #Live-demo with glob pattern
sel = run_dir.select('*/*DET/*', 'image.*')
sel.all_sources

Out[14]: frozenset({'SPB_DET_AGIPD1M-1/DET/0CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/10CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/11CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/12CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/13CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/14CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/15CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/1CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/2CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/3CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/4CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/5CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/6CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/7CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/8CHO:xtdf',
                    'SPB_DET_AGIPD1M-1/DET/9CHO:xtdf'})
```

Data can be access just like *previously* mentioned:

```
In [16]: #Live-demo (iterator example)
for train_id, train_data in sel.trains():
    #print(train_id, train_data.keys())
    break
train_id, train_data = sel.train_from_index(5)
train_id, train_data.keys()
```

```
Out[16]: (198425246,
dict_keys(['SPB_DET_AGIPD1M-1/DET/14CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/13CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/12CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/11CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/10CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/9CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/8CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/7CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/6CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/5CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/4CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/3CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/2CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/1CHO:xtdf', 'SPB_DET_AGIPD1M-1/DET/0CHO:xtdf']))
```

1.6 What if I don't want to fire up a python console just to quickly check something?

The command `lsxfel` is a tool that provides some basic functionality using only the command line.

```
$: lsxfel /gpfs/xfel/exp/XMPL/201750/p700000/raw/r0273/  
r0273 : Run directory
```

```
# of trains:      156  
Duration:         0:00:15.500000  
First train ID: 198425241  
Last train ID:   198425396
```

```
16 detector modules (SPB_DET_AGIPD1M-1)  
  e.g. module SPB_DET_AGIPD1M-1 0 : 512 x 128 pixels  
  176 frames per train, 27456 total frames
```

```
2 instrument sources (excluding detectors):  
  - SA1_XTD2_XGM/XGM/DOOCS:output  
  - SPB_XTD9_XGM/XGM/DOOCS:output
```

```
13 control sources:  
  - ACC_SYS_DOOCS/CTRL/BEAMCONDITIONS  
  - SA1_XTD2_XGM/XGM/DOOCS  
  - SPB_IRU_AGIPD1M/PSC/HV  
  - SPB_IRU_AGIPD1M/TSENS/H1_T_EXTHOUS  
  - SPB_IRU_AGIPD1M/TSENS/H2_T_EXTHOUS  
  - SPB_IRU_AGIPD1M/TSENS/Q1_T_BLOCK  
  - SPB_IRU_AGIPD1M/TSENS/Q2_T_BLOCK  
  - SPB_IRU_AGIPD1M/TSENS/Q3_T_BLOCK  
  - SPB_IRU_AGIPD1M/TSENS/Q4_T_BLOCK  
  - SPB_IRU_AGIPD1M1/CTRL/MC1  
  - SPB_IRU_AGIPD1M1/CTRL/MC2  
  - SPB_IRU_VAC/GAUGE/GAUGE_FR_6  
  - SPB_XTD9_XGM/XGM/DOOCS
```

1.7 If something goes wrong

A common source of errors is an invalid structure of the created data files. Hence a useful starting point to debug any errors is to check for valid file structures using the `karabo-data-validate` command:

```
$: karabo-data-validate /gpfs/xfel/exp/XMPL/201750/p700000/raw/r0273  
Checking run directory: /gpfs/xfel/exp/XMPL/201750/p700000/raw/r0273  
No problems found
```

The command checks if:

- All `.h5` files in a run can be opened, and the run contains at least one usable file.

- The list of train IDs in a file has no zeros except for padding at the end.
- Each train ID in a file is greater than the one before it.
- The indexes do not point to data beyond the end of a dataset.
- The indexes point to the start of the dataset, and then to successive chunks for successive trains, without gaps or overlaps between them.

1.8 Scenario II: Extracting 1D Data

Photon flux *time-series* (by trainID)


- How can 1D data be extracted and plotted?

The `get_series` method can extract a series across trainID's for a given device and property:

```
In [18]: #Live-demo
ph_flux = run_dir.get_series('SA1_XTD2_XGM/XGM/D00CS', 'pulseEnergy.photonFlux.value')
ph_flux.head()
```

```
Out[18]: trainId
198425241      500.519470
198425242      500.519470
198425243      502.727203
198425244      502.727203
198425245      504.070953
Name: SA1_XTD2_XGM/XGM/D00CS/pulseEnergy.photonFlux, dtype: float32
```

```
In [19]: #Live-demo
fig = plt.figure(figsize=(10,5))
_ = ph_flux.plot()
```



- Multiple devices and properties can be can be combined using the `get_dataframe` method:

```
In [20]: #Live-demo
fluxes_pos = run_dir.get_dataframe(fields=[("_XGM/*", "*.photonFlux"),
                                           ("*_XGM/*", "*.i[xy]Pos")])

fluxes_pos.head(10)
```

```
Out[20]:          SPB_XTD9_XGM/XGM/D00CS/pulseEnergy.photonFlux  \
trainId
198425241                404.392822
198425242                404.392822
198425243                404.392822
198425244                404.392822
198425245                404.392822
198425246                404.392822
198425247                404.665680
198425248                404.665680
198425249                406.852539
198425250                406.852539
```

```
          SPB_XTD9_XGM/XGM/D00CS/beamPosition.iyPos  \
trainId
198425241                -3.121433
198425242                -3.121433
198425243                -3.121433
198425244                -3.090523
198425245                -3.090523
198425246                -3.090523
198425247                -3.090523
198425248                -3.090523
198425249                -3.090523
198425250                -3.090523
```

```
          SPB_XTD9_XGM/XGM/D00CS/beamPosition.ixPos  \
trainId
198425241                5.512009
198425242                5.512009
198425243                5.512009
198425244                5.528512
198425245                5.528512
198425246                5.528512
198425247                5.528512
198425248                5.528512
198425249                5.528512
198425250                5.528512
```

```
          SA1_XTD2_XGM/XGM/D00CS/pulseEnergy.photonFlux  \
trainId
```

198425241	500.519470
198425242	500.519470
198425243	502.727203
198425244	502.727203
198425245	504.070953
198425246	504.070953
198425247	504.070953
198425248	504.070953
198425249	504.070953
198425250	505.071930

	SA1_XTD2_XGM/XGM/DOOCS/beamPosition.iyPos \
trainId	
198425241	0.315761
198425242	0.315761
198425243	0.315761
198425244	0.341187
198425245	0.341187
198425246	0.341187
198425247	0.341187
198425248	0.341187
198425249	0.341187
198425250	0.341187

	SA1_XTD2_XGM/XGM/DOOCS/beamPosition.ixPos
trainId	
198425241	1.293711
198425242	1.293711
198425243	1.293711
198425244	1.336566
198425245	1.336566
198425246	1.336566
198425247	1.336566
198425248	1.336566
198425249	1.336566
198425250	1.336566

1.8.1 What is *pandas*

get_dataframe and *get_series* return **pandas** objects which are extremely useful for extensive data analysis tasks. More information is available under <https://pandas.pydata.org>.

1.9 Scenario III: Getting data with multiple values per train

XGM intensity data is pulse resolved

- How can this 2D data be extracted and plotted?


the `get_array` method returns a data array that contains more than one value per train. XGM intensity data is pulse resolved and serves as an example.

```
In [21]: xgm_intensity = run_dir.get_array('SA1_XTD2_XGM/XGM/D00CS:output', 'data.intensityTD',  
                                           extra_dims=['pulseId'])
```

```
xgm_intensity
```

```
Out[21]: <xarray.DataArray (trainId: 156, pulseId: 1000)>  
array([[ 957.0532 , 1026.0005 ,  949.8755 , ...,    0.        ,    0.        ,  
         0.        ],  
       [ 763.8806 ,  794.2738 ,  868.2455 , ...,    0.        ,    0.        ,  
         0.        ],  
       [ 859.37    ,  995.1641 ,  838.5669 , ...,    0.        ,    0.        ,  
         0.        ],  
       ...,  
       [ 945.2731 ,  812.4336 ,  839.45654, ...,    0.        ,    0.        ,  
         0.        ],  
       [ 903.26855,  940.15125,  953.9436 , ...,    0.        ,    0.        ,  
         0.        ],  
       [ 944.08386,  949.549   ,  861.7509 , ...,    0.        ,    0.        ,  
         0.        ]], dtype=float32)  
Coordinates:  
  * trainId  (trainId) uint64 198425241 198425242 ... 198425395 198425396  
Dimensions without coordinates: pulseId
```

```
In [22]: #Live-demo  
fig = plt.figure(figsize=(10,8))  
_ = plt.imshow(xgm_intensity[:, :120].T,  
               extent=(xgm_intensity.trainId[0], xgm_intensity.trainId[-1], 1, 120),  
               origin='lower', cmap='RdYlBu_r')
```



a labeled array (**xarray**) is returned. More information on labeled arrays can be found on <http://xarray.pydata.org>

1.10 Scenario IV: Combine data and save new dataset to individual files

Suppose you want to compare data across different runs or combine various selections and then save the combined dataset to a new file.

This task could be realized with the **combine** and **write** method. This involves two different steps:

- combining
- file creation

```
In [20]: #Live-demo file combining
xgm_intensities = run_dir.select('*/XGM/*','data.intensityTD')
xgm_pos = run_dir.select('*/XGM/*', '*Pos')
xgm_union = xgm_pos.union(xgm_intensities)
```

Yet another convenient way of showing the content of a selection is the **selection** method:

```
In [21]: xgm_union.selection
```



```
Out[21]: {'SA1_XTD2_XGM/XGM/D00CS': {'beamPosition.ixPos.value',
    'beamPosition.iyPos.value'},
    'SA1_XTD2_XGM/XGM/D00CS:output': {'data.intensityTD'},
    'SPB_XTD9_XGM/XGM/D00CS': {'beamPosition.ixPos.value',
    'beamPosition.iyPos.value'},
    'SPB_XTD9_XGM/XGM/D00CS:output': {'data.intensityTD'}}
```

```
In [ ]: xgm_union.write('/gpfs/xfel/data/scratch/xmpl/xmpl_xgm_subset.hdf5')
```

```
$: file /gpfs/xfel/data/scratch/xmpl/xmpl_xgm_subset.hdf5
```

```
gpfs/xfel/data/scratch/xmpl/xmpl_xgm_subset.hdf5: Hierarchical Data Format (version 5) data
```

1.11 How to install/get karabo-data?

Karabo-data is available on GitHub and there are multiple ways to install it: * it is already available on the max-well cluster nothing has to be done * it can be install using *pip* &rarr `pip install (--user) karabo_data` (preferred way) * the latest version could be downloaded from GitHub &rarr `git clone https://github.com/European-XFEL/karabo_data.git`

```
In [ ]:
```